

# Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines  
and Vacuum Cold traps

## Lablyo Mini Freeze drier

### UK Manufactured

Bench top freeze drier with temperature display.  
20cm diameter x 25cm deep ice condenser  
chamber with 3kg ice capacity. 2kg ice collection  
in 24 hours.  
Condenser temperature  $-55^{\circ}\text{C}$ .  
Mobile temperature probe for product sampling.

### Part No's

Lablyo Mini	FIT/LYO/46/1
Vacuum pump for Lablyo	FIT/LYO/59/0

Can be used with the following freeze drying  
accessories:

### Part No's

Acrylic drying chamber	FIT/LYO/01/0
6 height shelf rack	FIT/LYO/03/0
6 heated shelf rack	FIT/LYO/19/0
4 shelf vial closing device	FIT/LYO/22/0
8 port flask manifold	FIT/LYO/13/0
48 port ampoule manifold	FIT/LYO/25/0
Spin Freezer	FIT/LYO/28/0



Apparatuur voor uw laboratorium  
De scherpste prijs voor de beste kwaliteit