

# Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines  
and Vacuum Cold traps

## Lablyo Plus Freezedrier

### UK Manufactured

This unique freeze drier is suitable for many applications in numerous disciplines. The **Lablyo Plus** offers standard features not normally found in bench top freeze driers. These include:

These include:

Stainless steel drying chamber,  
30cm diameter x 37cm deep with  
cooling and heating control, -40°C to +40°C

6 removable shelves inside the drying  
chamber in cake stand configuration.

Stainless steel condenser chamber,  
20cm diameter x 22cm deep. 3kg ice  
Capacity, 1kg ice collection per 24 hours  
-55°C condenser chamber temperature.

Electric de-frost system on condenser  
Chamber. 1.5hp none CFC refrigeration  
System. Condenser chamber top flange  
can accommodate other accessories  
from our range.

Touch screen with microprocessor control of temperatures and vacuum that are displayed on graph in real time. Mobile temperature probes in drying chamber for product sampling. The Lablyo Plus comes with clear acrylic lids for observing the process.

The separate condenser chamber and electric de-frost capability of the Lablyo Plus ensures maximum output and minimum turnaround time. The separate product and condenser chambers eliminate the need to pre-freeze samples and their temperature can be controlled during freeze-drying. Overall Dimensions: Height 57cm x length 92cm x depth x 52cm

The Lablyo Plus requires a vacuum pump and connection hose that is available separately.

### Part No's.

**Lablyo Plus**

**FIT/LYO/53/1**

**Vacuum pump**

**FIT/LYO/59/0**

**Vacuum hose**

**FIT/LYO/37/0**



### FROZEN IN TIME Ltd

York Road Industrial Park, Sheriff Hutton, York, YO60 6RZ England.

TEL. 01347 878158.

FAX. 01347 878303.

E-MAIL. [info@freezedriers.com](mailto:info@freezedriers.com)

WEBSITE. [www.freezedriers.com](http://www.freezedriers.com)