


DS-FS15E

desmonscientific.com/





BLOOD BANK CONTROLLER FEATURES

End user interface	Display LED
Resolution (°C)	0,1
Maximum Number of Probes	3
Internal Led light	yes
Acoustic alarm	yes
Hot gas advanced defrost	yes

ALARMS AND FAILURES (ACOUSTIC AND VISUAL)

HACCP function	yes
"Power failure" (back up battery) Alarm	yes
"High / low temperature" Alarm	yes
"Open door"Alarm	yes
"Sensor fault"Alarm	yes
Clogged condenser Alarm	yes
Temperature datalogger with USB port	yes

CONNECTIONS

RS 485 interface	optional
Porta Ethernet (LAN)	optional

SAFETY

Switch on/off protected by password	yes
-------------------------------------	-----

Dry Contact for alarms remoting	yes
---------------------------------	-----

Freezing and overheating safety device	yes
--	-----

OPTIONALS

Rechargeable back-up batteries	optional
--------------------------------	----------

STRUCTURE

Exterior dimensions WxDxH (mm)	600x600x1335
--------------------------------	--------------

Capacity (L)	150
--------------	-----

N of blood bags (450ml)	96
-------------------------	----

Insulation (mm)	60
-----------------	----

Temperature range (°C)	+4°C
------------------------	------

Max. ambient temperature (°C)	43
-------------------------------	----

Type of cooling unit	Ventilated
----------------------	------------

Power supply (V)	220V/240V
------------------	-----------

Frequency (Hz)	50
----------------	----

Power consumption (W)	245
-----------------------	-----

Weight (KG)	83
-------------	----

Shipping weight (KG)	94
----------------------	----

Shipping dimensions WxDxH (mm)	700x680x1720
--------------------------------	--------------

No. of compressors.	1
---------------------	---

Doors	1 glass door
-------	--------------

Automatic (smart) defrost and evaporation	yes
---	-----

Castors	yes
---------	-----

Chart recorder	yes
----------------	-----

Drawers (pcs)	3 s/steel drawers
Lock	yes
Interior material	s/steel AISI304
Exterior material	White powder coated steel
Gas	R404a

Price: € 0,00